IT - Cura Files SS e Ni-Ti K-File, H- File, K-Reamers

Strumenti canalari manuali

Descrizione

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Manual root canal instruments

Description

Manual root canal instruments made of steel are used in the probing and initial stages of root canal instrumentation. Preparation involves a widening of the root canal. Enlargement can be done by working in the crown-apical direction (crown-down technique) or in the apico-coronal direction (step-back technique) or both. Manual root canal instruments made of Nickel-Titanium have an S-shaped cross-section and allow rapid advancement in the canal with simultaneous removal of debris, without exerting pressure on the side walls of the canal. NI-Ti files are more flexible and can adapt to the curvature of the canal.

Indications

Only to be used in dental practices or in hospitals by qualified personnel.

Field of application

K-Files, H-Files, K-Reamers are used for glide path (No. 08 to 015) and endodontic shaping (all other diameters) during non-surgical endodontic treatment. K-Files 008 and 010 are also used in exploration, initial advancement and permeabilisation of the root canal.

Packaging

Non-sterile device in blister packs of 6 pieces (assorted or individual sizes).

Recommended procedure

Operating protocol K-File, H-File and K-Reamers:

- Locate the canal.
- Gently perform a passive probing of the coronal 2/3 with a hand instrument, with extreme gentleness and slight oscillation and determine the working length.
- If resistance is encountered, stop immediately. A file should never be forced.
- Before continuing, it is necessary to increase the coronal taper and probe the remaining section of the canal with a smaller hand file.
- The cutting action can be achieved either by filing (rasp action) or by a reaming. In a filing action, the tool is brought into the channel and then withdrawn, without rotation. With a reaming action, the file tends to get stuck in the dentine more easily and must therefore be handled with great care.
- Files can also be used in a push-pull action to widen certain canals with a filing action curved canals, as well as the ovoid section of wide canals.
- Irrigate copiously and clean the instruments constantly to prevent debris from being compacted at the foramen or in the periapical tissues.
- Then proceed with the most suitable operating sequence for the case to be treated.

During the endodontic sequence, irrigate repeatedly with hypochlorite and E.D.T.A. in order to remove all necrotic tissue and the organic component of the dentinal sludge.

Storage

Store at room temperature, away from moisture and direct light.

Precautions

- Before use, it is recommended to check with patients for any allergy to raw materials.
- It is recommended to always use surgical gloves for personal protection against bacterial contamination.
- Failure to do so may result in cross-infection.

Warnings

- Single-use device, re-use is not recommended to limit problems with cross-infection and/or breakage.
- Do not use the blister if opened or damaged.
- Never force the instrument into the canal.
- Do not use after the indicated expiry date.
- You must report any serious incidents occurring in relation to the medical device you have supplied to the manufacturer and to the competent authority in your Member State.

Contraindications

Do not use in patients with known allergic sensitivity to raw materials.

Symbols

Product code

Lot number

Expiry date

Do not reuse

Keep dry

Keep away from sunlight

Non-sterile

IFU

Packing unit

Do not use if packaging is damaged

Medical device

Unique Device Identification

Stainless Steel

Nickel-Titanium

Manufacturer

Warranty

Dental World shall replace the product if it proves to be defective.

Dental World does not accept liability for any damage or loss, direct or consequential, resulting from careless use or inability to use the product described. Before use, it is the user's responsibility to determine the appropriateness of the product for its use. The user assumes responsibility for careless use.

Dental World S.r.l.
Via Antichi Pastifici, 15 - 70056 Molfetta (BA)
Tel. +39 080 3381004 - Fax. +39 080 3386617
e-mail: info@dental-world.it / web: www.dental-world.it

Curaden AG - Switzerland
Amlehnstrasse 22 - 6010 Kriens
www.curaden.com / info@curaden.ch

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