# IT - Cura Pat

## Strumenti canalari meccanici

### Descrizione

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Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

Curaden AG - Switzerland

Amlehstrasse 22 - 6010 Kriens

www.curaden.com / info@curaden.ch

# EN - Cura Pat

## Mechanical root canal instruments

### Description

Cura Pat root canal instruments are designed for greater flexibility and resistance to cyclical fatigue. They feature a gradual increase in diameter and a rounded tip that avoids steps and zips and optimises the distance between blades to increase their strength and at the same time their ability to remove debris. Also suitable in calcified and long canals. These new rotary instruments are designed to quickly and safely create mechanical pre-flaring and glide-path, thus eliminating the last manual step.

### Indications

Only to be used in dental practices or in hospitals by qualified personnel.

### Field of application

Shaping and cleaning of root canals.

### Measures

Cura Pat								
	Ø	%	N/cm	RPM	21 mm	25 mm	31 mm	
● 013	13	2	3	300	155900460	155900470	155900480	
○ 016	16	2	3	300	155900461	155900471	155900481	
● 019	19	2	3	300	155900462	155900472	155900482	
013-019	--	--	3	300	155900463	155900473	155900483	

### Packaging

Non-sterile device in blister packs of 6 pieces (assorted or individual sizes).

### Recommended procedure

#### Operating protocol:

- In the absence of coronal interference, probe with an ISO 010 steel K-File until it moves freely in the canal (glidepath) and thus verify the working length.
- After establishing the patency of the apical foramen and confirming the working length, use the 3 Cura Pat files at the same length, thus performing pre-flaring in a few seconds.
- If the working length can be taken with a K-File ISO 015, go directly to the Cura Pat 0.16 mm.
- If the working length can be taken with a K-File ISO 020, one goes directly to Cura Pat 0.19 mm.
- After obtaining a 19-gauge foramen, we can safely take the rotary system we are accustomed to using, such as Cura Shaper or others, to the same working length and safely perform the correct shaping.
- After the use of each instrument, it is advisable to perform a thorough irrigation with sodium hypochlorite and viscous chelator as a lubricant.
- In the event that the canal cannot be probed up to the foramen from the beginning, due to the presence of coronal If the canal cannot be probed all the way to the foramen from the beginning, due to coronal interference or very sharp curves in the apical third, the Cura Pat is used in step back, taking it to where the canal accommodates it without engaging the tip, and then probing is performed with the K-File 10 pre-curved, take the working length and proceed as described above.

Throughout the endodontic sequence, repeatedly irrigate with hypochlorite and E.D.T.A. in order to remove all necrotic tissue and the organic component of the dentinal slurry.

### Storage

Store at room temperature, away from moisture and direct light.

### Precautions

Before use, it is recommended to check with patients for any allergy to raw materials.

It is recommended to always use surgical gloves for personal protection against bacterial contamination. Failure to do so may result in cross-infection.


### Warnings

- Single-use device, re-use is not recommended to limit problems with cross-infection and/or breakage.
- Do not use the blister if opened or damaged.
- Never force the instrument into the canal.
- Do not use after the indicated expiry date.
- You must report any serious incidents occurring in relation to the medical device we have supplied to the manufacturer and to the competent authority in your Member State.

### Contraindications

Do not use in patients with known allergic sensitivity to raw materials.

### Symbols

### Warranty

Dental World shall replace the product if it proves to be defective.

Dental World does not accept liability for any damage or loss, direct or consequential, resulting from careless use or inability to use the product described. Before use, it is the user's responsibility to determine the appropriateness of the product for its use. The user assumes responsibility for use careless use.

Dental World S.r.l.  
Via Antichi Pastifici, 15 - 70056 Molfetta (BA)  
Tel. +39 080 3381004 - Fax. +39 080 3386617  
e-mail: info@dental-world.it / web: www.dental-world.it

Curaden AG - Switzerland  
Amlehnstrasse 22 - 6010 Kriens  
www.curaden.com / info@curaden.ch

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