IT - ELI Flex Evolution

Strumenti canalari meccanici

Descrizione

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EN - ELI Flex Evolution

Mechanical root canal instruments

Description

ELI Flex Evolution root canal instruments are thermally activated Ni-Ti and Coated files designed for complete root canal preparation. The kit includes an opening instrument, the ELI Pilot N.H.A. glide-path instrument and four rotary files for shaping and finishing all types of canals.

Indications

For use only in a dental practice or hospital environment by qualified personnel.

Field of application

Glide path, shaping and cleaning root canals.

Measures

ELI Flex Evolution							
	Ø	%	N/cm	RPM	21 mm	25 mm	31 mm
R	20	10	2,5	350	17 mm REF: 7301028		
● D1	20	4	2,5	150-350	7301029	7301039	7301049
● D2	25	4	2,5	150-350	7301030	7301040	7301050
● D3	30	4	2,5	150-350	7301031	7301041	7301051
● D4	25	6	2,5	150-350	7301032	7301042	7301052
R-D4	--	--	2,5	--	7301033	7301043	7301053
● D5	35	4	2,5	150-350	7301034	7301044	7301054
● D6	40	4	2,5	150-350	7301035	7301045	7301055
● D7	50	4	2,5	150-350	7301036	7301046	7301056
D5-D7	--	--	2,5	--	7301037	7301047	7301057

Packaging

Non-sterile device in blister packs of 6 pieces (assorted or individual sizes).

Recommended procedure

Operating protocol:

- Locate the canal.
- Gently perform a passive probing of the coronal 2/3 with a hand instrument ISO 010 or 015 (mechanical files such as ELI Pat can be used alternatively).
- Using the R instrument, create access to the canal orifice and irrigate.
- Use sodium hypochlorite and viscous chelator as a lubricant to irrigate the canal.
- To achieve a reproducible glide-path before enlargement and to reduce the workload of the shaping instruments and the risk of breakage, work gently with the ELI Pilot N.H.A. rotary instrument (ISO 016).
- Shape with D1 until working length is reached, irrigate.
- Use D2 to reconfirm the working length, recapitulate and irrigate again, especially in more curved canals.
- For large canals use D3/D4 without oscillating movement.
- For finishing rather large root canals use instruments D5, D6 and D7.

Throughout the endodontic sequence, repeatedly irrigate with hypochlorite and E.D.T.A. in order to remove all necrotic tissue and the organic component of the dentinal slurry.

Storage

Store at room temperature, away from moisture and direct light.

Precautions

Before use, it is recommended to check with patients for any allergy to raw materials.

It is recommended to always use surgical gloves for personal protection against bacterial contamination. Failure to do so may result in cross-infection.


Warnings

- Single-use device, re-use is not recommended to limit problems with cross-infection and/or breakage.
- Do not use the blister if opened or damaged.
- Never force the instrument into the canal.
- Do not use after the indicated expiry date.
- You must report any serious incidents occurring in relation to the medical device we have supplied to the manufacturer and to the competent authority in your Member State.

Contraindications

Do not use in patients with known allergic sensitivity to raw materials.

Symbols

Warranty

Dental World shall replace the product if it proves to be defective.

Dental World does not accept liability for any damage or loss, direct or consequential, resulting from careless use or inability to use the product described. Before use, it is the user's responsibility to determine the appropriateness of the product for its use. The user assumes responsibility for use careless use.

Dental World S.r.l.
Via Antichi Pastifici, 15 - 70056 Molfetta (BA)
Tel. +39 080 3381004 - Fax. +39 080 3386617
e-mail: info@dental-world.it / web: www.dental-world.it

Eli-Dent Group S.p.A.
Via Guido Rossa, 24 A/B/C - 50032
Borgo San Lorenzo (FI) - www.elidentgroup.it
Tel. +39 055 600545 - Fax +39 055 601069

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