Dabliu Files SS e Ni-Ti

K-File, H- File, K-Reamers

Strumenti canalari manuali

Descrizione

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Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

Dabliu Files SS e Ni-Ti

K-File, H- File, K-Reamers

Manual canal instruments

Description

Manual stainless steel canal instruments are used in sampling phases and in the initial stages of the canal instrumentation. Preparation involves an enlargement of the root canal. Enlargement can take place by working in crown-apical direction (crown-down technique) or apical-coronal direction (step-back technique) or both.

The manual canal instruments in Nickel-Titanium, present S-section and allow a quick advance in the canal with simultaneous removal of debris, without exerting pressure on the side walls of the canal. Ni-Ti Files are more flexible and able to adapt to canal curvatures.

Indications

To be used only in dental clinics or hospitals.

Application field

K-Files, H-Files, K-Reamers tools are used for the glide path (numbers 08 to 015) and for the endodontics shaping (all other diameters) during non-surgical endodontic treatments. K-Files 008 and 010 are used also for exploration, in the initial advancement and in the permeability of the root canal.

Packaging

Blister 6 pcs (assorted or individual sizes).

Recommended procedure

Operative protocol K-File, H-File and K-Reamers:

- Identify the canal.
- Gently carry out a passive sounding of 2/3 coronals with a manual file, with extreme delicacy and a slight swing and determine the working length.
- If resistance is encountered, you must stop immediately. A file should never be forced.
- Before proceeding, you need to increase the coronal taper and probe the remaining stretch of the canal with a smaller manual file.
- The cutting action can be achieved either by a filing (scraping action), either with an action of boring. In a filing motion, the instrument is brought into the canal and then withdrawn, without rotation. With a boring action the file tends to get stuck in the dentine more easily and should therefore be handled with great caution.
- The files can also be used in push-drive (traction) to widen by filing certain curved canals, as well as the

avoid section of wide canals.

- Thoroughly irrigate and constantly clean the instrument to prevent debris from being compacted in the foramen or in the periapical tissues.
- Then proceed with the operational sequence more appropriate to the case to be dealt with.

Throughout the endodontic sequence, irrigate repeatedly with hypochlorite (Endo Clor) and E. D. T. A. (Canal Prep) in order to remove all necrotic tissue and organic component of the smear layer.


Storage

Store at room temperature, away from direct light sources.

Cautions

- The instruments are recommended for single use.
- If used repeatedly they may not be efficient and/or cause inconvenience.
- Do not use the blister if opened or damaged.

Symbols

Warranty

Dental World will replace product that is proven to be defective. Dental World does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risks and liability in connection therewith.

 Dental World S.r.l.

Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it