Pro Pat

Strumenti canalari meccanici

Descrizione

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Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

0476

Pro Pat

Mechanical canal instruments

Description

Pro Pat instruments are designed to achieve greater flexibility and resistance to cyclical fatigue, characterized by a gradual increase in diameter and a rounded tip that avoids steps, zip and optimizes the distance between the blades, to increase their strength and their ability to remove debris. Also suitable in calcified and long canals. These new rotating instruments are designed to create quickly and safely pre-flaring and mechanical glide-path, thus eliminating the last manual step.

Indications

To be used only in a dental clinics or in hospitals.

Application field

Shaping and cleaning of root canal.

Measures

PRO PAT							
	Ø	%	N/cm	RPM	21 mm	25 mm	31 mm
● 013	13	2	3	300	144900460	144900470	144900480
○ 016	16	2	3	300	144900461	144900471	144900481
● 019	19	2	3	300	144900462	144900472	144900482
013-019	--	--	3	300	144900463	144900473	144900483

Packaging

6 pcs (assorted or in individual sizes).

Recommended procedure

Operative protocol:

- In the absence of coronal interference, perform a sounding with a K-File steel file ISO 010 until it moves freely in the canal and check the working length.
- After establishing the pervity of the apical foramen and confirming the working lenght, use the 3 Pro Pat files with the same length, executing in a few seconds the pre-flaring.
- If the working lenght can be taken with a K-File ISO 015, use directly to the Pro Pat 16 mm.
- If the l. l. can be taken with a K-File ISO 020, go directly to the Pro Pat 19 mm.
- After obtaining a 19-gauge foramen, you can use the same working lenght for the rotating system you usually use and safely execute the correct shaping.
- After using each instrument, it is advisable to irrigate abundantly with sodium hypochlorite (Endo Clor) and viscous chelator as a lubricant (Canal Prep).
- In the event that the canal is not probeable until the foramen from the beginning, due to coronal interference

or very marked curves in the third apical, Pro Pat is used in step back taking them to where the canal receives them without compromising the tip of the canal, and then the survey is carried out with the file K-File 10 pre-curved, it takes the working length and proceed as described above.

Throughout the endodontic sequence, irrigate repeatedly with hypochlorite (Endo Clor) and E. D. T. A. (Canal Prep) in order to remove all necrotic tissue and organic component of the smear layer.


Storage

Store at room temperature, away from direct light sources.

Preventive measures

- The instruments are recommended for single use.
- If used repeatedly they may not be efficient and/or cause inconvenience.
- Do not use the blister if opened or damaged.

Symbols

Warranty

Dental World will replace product that is proven to be defective. Dental World does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risks and liability in connection therewith. e la responsabilità per uso incauto.

Dental World S.r.l.

Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it