Pro Race H.A.

Strumenti canalari meccanici

Descrizione

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Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

Pro Race H.A.

Mechanical canal instruments

Description

Pro Race H.A. instruments are made of Ni-Ti C New Wire, which confers greater elasticity and higher strength than normal Ni-Ti instruments.

Indications

To be used only in dental clinics or in hospitals.

Application field

Shaping and cleaning of root canal.

Measures

PRO RACE <i>Heat Activation</i>					
	∅	%	21 mm	25 mm	31 mm
● R25	25	0,8	144900710	144900720	144900725
● R40	40	0,6	144900711	144900721	144900726
● R50	50	0,5	144900712	144900722	144900727
R25-R50	--	--	144900713	144900723	144900728

Packaging

Blister 6 pcs (assorted or individual sizes).

Recommended procedure

Operative protocol:

- Arrange an access to the root canal.
- Select the Pro Race instrument to use for the preparation of the canal.
- If the canal is considered narrow use the tool R25. If the channel is considered medium use the instrument R40. If the canal is considered wide use the R50 instrument.
- Preparation with the R25 tool: estimate the length of the root canal with the help of X-ray and adjust the stop on the instrument to about 2/3 of the estimated length.
- Preparation with the R40-R50 instrument: set the working length in electronic mode and adjust the stop on the instrument at that length.
- Irrigate the canal with E. D. T. A. (Canal Prep) and hypochlorite (Endo Clor).
- Insert the tool into the canal and move it slowly in the apical direction, applying a very light pressure, so the instrument can move easily. The range of the up and down movements (passages) shall not exceed 3 mm.
- After 3 steps, clean the instrument, irrigate the canal and check its pervity.
- Recapitulate until you reach 2/3 of the estimated working length.
- If the R25 instrument is used, electronically set working length. If you are using the R40 and the R50 instruments, after preparing 2/3 of the length of the canal, electronically check the established working

length.

- Next, recapitulate until you reach the entire working length.

Determine the apical diameter by inserting one ISO manual instrument superior to that of the Pro Race instrument working length. If the manual instrument is inserted up to 0.7 mm before the working length, no further extension of the apical area is required. If the manual instrument can be inserted up to working length, it is usually necessary to finish the preparation with a larger Pro Race instrument.

Storage

Store at room temperature, away from direct light sources.

Preventive measures

- The instruments are recommended for single use.
- If used repeatedly they may not be efficient and/or cause inconvenience.
- Do not use the blister if opened or damaged.
-


Storage

Store at room temperature, away from direct light sources.

Preventive measures

- The instruments are recommended for single use.
- If used repeatedly they may not be efficient and/or cause inconvenience.
- Do not use the blister if opened or damaged.

Symbols

Warranty

Dental World will replace product that is proven to be defective. Dental World does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risks and liability in connection therewith.

 Dental World S.r.l.

Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

 0425