IT - Pro Pat

Strumenti canalari meccanici

Descrizione

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Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

EN - Pro Pat

Mechanical root canal instruments

Description

Pro Pat root canal instruments are designed for greater flexibility and resistance to cyclical fatigue. They feature a gradual increase in diameter and a rounded tip that avoids steps and zips and optimises the distance between blades to increase their strength and at the same time their ability to remove debris. Also suitable in calcified and long canals. These new rotary instruments are designed to quickly and safely create mechanical pre-flaring and glide-path, thus eliminating the last manual step.

Indications

Only to be used in dental practices or in hospitals by qualified personnel.

Field of application

Shaping and cleaning of root canals.

Measures

PRO PAT							
	Ø	%	N/cm	RPM	21 mm	25 mm	31 mm
● 013	13	2	3	300	144900460	144900470	144900480
○ 016	16	2	3	300	144900461	144900471	144900481
● 019	19	2	3	300	144900462	144900472	144900482
013-019	--	--	3	300	144900463	144900473	144900483

Packaging

Non-sterile device in blister packs of 6 pieces (assorted or individual sizes).

Recommended procedure

Operating protocol:

- In the absence of coronal interference, probe with an ISO 010 steel K-File until it moves freely in the canal (glidepath) and thus verify the working length.
- After establishing the patency of the apical foramen and confirming the working length, use the 3 Pro Pat files at the same length, thus performing pre-flaring in a few seconds.
- If the working length can be taken with a K-File ISO 015, go directly to the Pro Pat 0.16 mm.
- If the working length can be taken with a K-File ISO 020, one goes directly to Pro Pat 0.19 mm.
- After obtaining a 19-gauge foramen, we can safely take the rotary system we are accustomed to using, such as Pro Shaper or others, to the same working length and safely perform the correct shaping.
- After the use of each instrument, it is advisable to perform a thorough irrigation with sodium hypochlorite (Endo Clor) and viscous chelator as a lubricant (Canal Prep).
- In the event that the canal cannot be probed up to the foramen from the beginning, due to the presence of coronal If the canal cannot be probed all the way to the foramen from the beginning, due to coronal interference or very sharp curves in the apical third, the Pro Pat is used in step back, taking it to where the canal accommodates it without engaging the tip, and then probing is performed with the K-File

10 pre-curved, take the working length and proceed as described above.

Throughout the endodontic sequence, repeatedly irrigate with hypochlorite (Endo Clor) and E.D.T.A. (Canal Prep) in order to remove all necrotic tissue and the organic component of the dentinal slurry.

Storage

Store at room temperature, away from moisture and direct light.

Precautions

Before use, it is recommended to check with patients for any allergy to raw materials. It is recommended to always use surgical gloves for personal protection against bacterial contamination. Failure to do so may result in cross-infection.


Warnings

- Single-use device, re-use is not recommended to limit problems with cross-infection and/or breakage.
- Do not use the blister if opened or damaged.
- Never force the instrument into the canal.
- Do not use after the indicated expiry date.
- You must report any serious incidents occurring in relation to the medical device we have supplied to the manufacturer and to the competent authority in your Member State.

Contraindications

Do not use in patients with known allergic sensitivity to raw materials.

Symbols

Warranty

Dental World shall replace the product if it proves to be defective.

Dental World does not accept liability for any damage or loss, direct or consequential, resulting from careless use or inability to use the product described. Before use, it is the user's responsibility to determine the appropriateness of the product for its use. The user assumes responsibility for use careless use.

 Dental World S.r.l.
Via Antichi Pastifici, 15 - 70056 Molfetta (BA)
Tel. +39 080 3381004 - Fax. +39 080 3386617
e-mail: info@dental-world.it / web: www.dental-world.it

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