IT - Dabliu Files SS e Ni-Ti K-File, H- File, K-Reamers

Strumenti canalari manuali

Descrizione

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Via Antichi Pastifici, 15 - 70056 Molfetta (BA)

Tel. +39 080 3381004 - Fax. +39 080 3386617

e-mail: info@dental-world.it / web: www.dental-world.it

EN - Dabliu SS and Ni-Ti Files

K-File, H- File, K-Reamers

Manual root canal instruments

Description

Manual root canal instruments made of steel are used in the probing and initial stages of root canal instrumentation. Preparation involves a widening of the root canal. Enlargement can be done by working in the crown-apical direction (crown-down technique) or in the apico-coronal direction (step-back technique) or both. Manual root canal instruments made of Nickel-Titanium have an S-shaped cross-section and allow rapid advancement in the canal with simultaneous removal of debris, without exerting pressure on the side walls of the canal. NI-TI files are more flexible and can adapt to the curvature of the canal.

Indications

Only to be used in dental practices or in hospitals by qualified personnel.

Field of application

K-Files, H-Files, K-Reamers are used for glide path (No. 08 to 015) and endodontic shaping (all other diameters) during non-surgical endodontic treatment. K-Files 008 and 010 are also used in exploration, initial advancement and permeabilisation of the root canal.

Packaging

Sterile device in blister packs of 6 pieces (assorted or individual sizes).

Recommended procedure

Operating protocol K-File, H-File and K-Reamers:

- Locate the canal.
- Gently perform a passive probing of the coronal 2/3 with a hand instrument, with extreme gentleness and slight oscillation and determine the working length.
- If resistance is encountered, stop immediately. A file should never be forced.
- Before continuing, it is necessary to increase the coronal taper and probe the remaining section of the canal with a smaller hand file.
- The cutting action can be achieved either by filing (rasp action) or by a reaming. In a filing action, the tool is brought into the channel and then withdrawn, without rotation. With a reaming action, the file tends to get stuck in the dentine more easily and must therefore be handled with great care.
- Files can also be used in a push-pull action to widen certain canals with a filing action curved canals, as well as the ovoid section of wide canals.
- Irrigate copiously and clean the instruments constantly to prevent debris from being compacted at the foramen or in the periapical tissues.
- Then proceed with the most suitable operating sequence for the case to be treated.

During the endodontic sequence, irrigate repeatedly with hypochlorite (Endo Clor) and E.D.T.A. (Canal Prep) in order to remove all necrotic tissue and the organic component of the dentinal sludge.

Storage

The blister should only be opened at the time of surgery and under controlled aseptic conditions.

Store at room temperature, away from moisture and direct light.

Precautions

- Before use, it is recommended to check with patients for any allergy to raw materials.
- It is recommended to always use surgical gloves for personal protection against bacterial contamination.
- Failure to do so may result in cross-infection.


Warnings

- Single-use device, re-use is not recommended to limit problems with cross-infection and/or breakage.
- Do not use the blister if opened or damaged.
- If the package shows signs of tampering or damage due to transport or any other nature, the instruments contained therein may not be used as they may have lost sterility.
- Never force the instrument into the canal.
- Do not use after the indicated expiry date.
- You must report any serious incidents occurring in relation to the medical device you have supplied to the manufacturer and to the competent authority in your Member State.

Contraindications

Do not use in patients with known allergic sensitivity to raw materials.

Symbols

Warranty

Dental World shall replace the product if it proves to be defective.

Dental World does not accept liability for any damage or loss, direct or consequential, resulting from careless use or inability to use the product described. Before use, it is the user's responsibility to determine the appropriateness of the product for its use. The user assumes responsibility for careless use.

 Dental World S.r.l.
Via Antichi Pastifici, 15 - 70056 Molfetta (BA)
Tel. +39 080 3381004 - Fax. +39 080 3386617
e-mail: info@dental-world.it / web: www.dental-world.it

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