


PRO FLEX *N.H.A*


NANO-COATED HEAT ACTIVATION



PRO FLEX NHA

Pro Flex NHA is a systematics designed with a triangular cross section and a proprietary heat treatment that gives it excellent memory control and flexibility for bends up to 90°, adapting to even the most complex root canal anatomies.


The technology used in **Pro Flex NHA's** heat activation treatment and nanoparticles coating provide the instruments with **70% more flexibility** and up to **3 times more cyclic fatigue** resistance than any heat-treated instrument.




CONTROL
MEMORY



TRANGULAR CROSS SECTION
PRECISE AND SHARP BLADE




Flexural strength



Number of revolutions to breakage

Cyclic torsion resistance



Torque value at breakage

LEGEND :

- Pro Flex NHA
- Competitor 1
- Competitor 2



BEFORE


Radiographic evidence of root canal pathology.



AFTER

Radiograph of the perfectly treated and sealed root canal thanks to the precision of **Pro Flex NHA** instruments.

PRO FLEX N.H.A SYSTEMATICS



Coming Soon

The **Pro Flex NHA** systematics marks a new era in root canal treatment, distinguished by its ease of use. This systematics, makes treatments significantly easier, respecting the anatomy of the root canal by avoiding protrusions, weakening or perforations with unparalleled precision.

Created and developed exclusively by Dental World, **Pro Flex NHA** instruments offer professionals an intuitive and reliable alternative to traditional root canal instruments while ensuring restorative treatment excellence.

PACKAGING



Pro Flex NHA root canal instruments are packaged in 6-pieces blister packs, providing convenience, integrity and protection from contamination.

PRO FLEX NHA									
	∅	%	N/cm	RPM	21 mm	25 mm	31 mm	Cross-section	Match gutta point
R	20	10	2,5	350	17 mm REF: 144900550				-
○ 016	15	2-6	2,5	300	144900551	144900561	144900571		-
● D1	20	4	2,5	150-350	144900552	144900562	144900572		20/04
● D2	25	4	2,5	150-350	144900553	144900563	144900573		25/04
● D3	30	4	2,5	150-350	144900554	144900564	144900574		30/04
● D4	25	6	2,5	150-350	144900555	144900565	144900575		25/06
R-D4	--	--	2,5	--	144900556	144900566	144900576		-
● D5	35	4	2,5	150-350	144900557	144900567	144900577		35/04
● D6	40	4	2,5	150-350	144900558	144900568	144900578		40/04
● D7	50	4	2,5	150-350	144900559	144900569	144900579		50/04
D5-D7	--	--	2,5	--	144900560	144900570	144900580		-

Coming SOON

* Speed and torque are approximate and may vary depending on the device used and operator preference.
Legend: ∅: Diameter; %: Taper; N/cm: Torque; RPM: Speed



Manufactured by:

 **DENTAL WORLD S.r.l.**

Via Antichi Pastifici, 15 Molfetta (BA) - ITALY
Tel. +39 080 3381004 / Fax. +39 080 3386617
www.dental-world.it / info@dental-world.it