



# PROALGIN

Alginates for dental impressions



MADE IN ITALY

# PROALGIN ACCURATE


**PROALGIN Accurate** is the new generation alginate with **reproduction value up to 15µ**, used in any kind of operation where the **definition of detail** is an essential condition for the correct success of a prosthesis. It is also particularly suitable for fixed prostheses as it retains all the points of strength while not omitting a long dimensional stability.

## Characteristics


- Class A, type 1
- Easily mixable
- Thixotropic
- High elasticity
- High detail definition: 15 µm
- High dimensional stability (168 hours)
- Compatible with impression disinfectants
- Dust-free
- Mint flavor

## Applications

- Removable prostheses
- Antagonists in removable and fixed prostheses
- Dental provisionals
- Preliminary models
- Skeleton prostheses
- Positional impressions



## Working phases



PROALGIN ACCURATE - 450g bag product code: 001033N

# PROALGIN CHROM

**PROALGIN Chrom** is the first alginate that, with a **chromatic phase indicator**, revolutionised the world of dentistry: impressions with **high dimensional stability**, easy control of the working steps and quick setting time. Temperature variations that cause changes in processing times are automatically indicated in the colour shadings.

## Characteristics


- Class A, type 1
- Easily mixable
- Thixotropic
- High detail definition: 20 µm
- High dimensional stability (168 hours)
- Compatible with impression disinfectants
- Dust-free
- Mint flavor

## Applications

- Removable prostheses
- Antagonists in removable and fixed prostheses
- Preliminary models
- Study models
- Skeleton prostheses



## Working phases



PROALGIN CHROM - 450g bag product code: 001025N