

# RVG POSITIONER

Sistema di centratura per sensori digitali  
*Digital sensor's centering system*



## Descrizione

RVG Positioner è un prodotto realizzato in materiale plastico ed indicato per il posizionamento di pellicole per sensori digitali a raggi-X, al fine di ottenere una maggiore accuratezza nella procedura di ripresa radiografia.

## Indicazioni:

Utilizzato per fotografare i denti, aiutando il medico a trovare la posizione esatta per poter scattare delle foto. I pezzi si differenziano attraverso dei colori codificati, sono leggeri, semplici da usare e facili da montare. Tutti i supporti sono autoclavabili a 121°C.

## Description

RVG Positioner is a product made of plastic and suitable for positioning films for digital X-ray sensors in order to achieve greater accuracy in the X-ray imaging procedure.

## Indications:

Used for photographing teeth, helping the dentist find the exact position for taking pictures. The pieces are differentiated through color coded, are lightweight, simple to use and easy to assemble. The products are autoclavable at 121°C.



## Dati tecnici


### Technical data

## SET DA 3 PEZZI

- **Verde:** per premolare/molare superiore mascellare destro e mandibolare sinistro.
- **Blu:** per premolare mascellare sinistro e premolare/molare mandibolare destro.
- **Giallo:** per incisivo mascellare/mandibolare denti anteriori.

## 3 PIECES SET

- **Green:** for right maxillary and left mandibular upper premolar/molar.
- **Blue:** for left maxillary premolar and right mandibular premolar/molar.
- **Yellow:** for maxillary/mandibular incisor anterior teeth.



DENTAL WORLD S.r.l.

Via Antichi Pastifici, 15 - 70056 Molfetta (BA) ITALY  
tel. +39 080.33810.04 / fax. +39 080.338.6617  
[www.dental-world.it](http://www.dental-world.it) / [info@dental-world.it](mailto:info@dental-world.it)

PACKAGING - RVG POSITIONER REF: 00520040

- Confezione da 3 pezzi
- Istruzioni
- Prezzo listino: € 56,00 + IVA
- Box of 3 pcs
- IFU
- Price list: € 56,00 + IVA